

Truly interactive hearing



Get to know Yuu™

What makes Yuu™? In a word, interaction. Your clients experience exceptional listening in diverse environments and situations. They have the control to make real-time, personal refinements to that experience. And with each refinement, they teach Yuu their listening preferences for every situation.



So much more than your average remote control: Yuu is also available with Smart Control, offering added functionality and convenience.





SIGNATURE FEATURES

AutoPro4™	Experience superior automatic with faster, smoother transitions
Comfort–Clarity Balance	Industry first: enables clients to personalize levels of speech enhancement and noise reduction either on the hearing instrument or with Smart Control
Self learning	Gradually learns client preferences in the automatic program
Feedback Manager	Breakthrough phase canceller delivers more useable gain
Data Logging	Enhanced: provides data from Comfort-Clarity Balance and volume control in the automatic and manual programs
Easy DAI	Automatically senses and amplifies external audio devices on 13 BTE models
AntiShock™	Eliminates discomfort from sudden impulse noises

SMART CONTROL WITH LEARN NOW™

Smart Control	Industry first: remote control with Comfort-Clarity Balance, learnNow, program changes, volume control, and binaural control
LearnNow™	Industry first: instantly captures client preferences in the automatic program

ADDITIONAL FEATURES

Channels	20
Programs (Auto, Manual)	Automatic plus up to 3 manual
Processing Schemes	WDRC
Directional System	20 band adaptive directional
Speech Enhancement	Level dependent
Noise Reduction	20 bands
MyMusic™	Enhances the experience of listening to music by bringing out rich, full tones
Automatic telephone solution	Easy-t (custom and BTE)
Models/shell styles available	Custom CIC thru Full Shell Power Micro BTE 312 thru 13 HP
Micro BTE (312)	Moda™ II, Moxi™
Slim tube options	Moda II; 13 Standard BTE
Canal Receiver Technology (CRT)	Moxi xS receiver Moxi xP receiver
OnBoard™ Control	Choice of volume control or program button
Wind Noise Manager	✓
Insitu Audiogram	✓

Get to know Yuu
Visit www.unitronhearing.com

